

# Technical Data Sheet

## Application Sheet for Citifluor AF3 Mountant Solution

#17972-25, 17972-100

### Usage

AF3 is a mountant solution composed of phosphate buffered saline and an antifadent. It was specifically designed to stop the photobleaching of the fluorescein moiety of FITC labelled biological specimens. Its application is not, however, limited to FITC labeled materials and has been used with advantage with many other fluorochromes including rhodamines, Alexa dyes and DAPI. It is useful for examining tissue live cells although it should be noted that cell lysis will take place over a period of time. If this is a problem, use our AFR3 solution.

The solution has a pH of ~10. The solution should be pipetted onto the specimen and then a cover slip applied. Specimens mounted in AF3 solutions have been kept in this way for many months without suffering damage.

### Properties and Storage of AF3 Mountant Solution

The solutions have a water-white in appearance. They may be stored at room temperature and ideally between 5° and 15°C and out of strong sunlight. The cap of the bottle or if using the pipette supplied with the material, the cap which covers the pipette delivery point, should always be replaced after use as a matter of good practice. Samples stored under these conditions for 6 months have shown no apparent deterioration. If the AF3 solution is being used in an assay, a control experiment should always be carried out.